Deposition Checklist

Check sample sequence

Refer to UniProt

Check ligands

Are these correctly assigned?

✓ Prepare data for deposition

Generate PDBx/mmCIF-formatted data files to provide the most complete information about your structure

✓ Prepare data for deposition

Generate PDBx/mmCIF-formatted data files to provide the most complete


Output

Description

Validate your data
Use wwPDB's standalone
validation service

Deposit your structure using OneDep System

When you are ready



Online Resources

PDB_extract: pdb-extract.wwpdb.org
Validation: validate.wwpdb.org
Deposition: deposit.wwpdb.org

Further Information

Tutorials and FAQs for validation and deposition can be found at:

www.wwpdb.org/validation/validation-reports wwpdb.org/deposition/tutorial wwpdb.org/deposition/faq

Contact Information

validation@mail.wwpdb.org deposit-help@mail.wwpdb.org

Citation Information

Cite wwPDB:

Nature Structural Biology **10**, 980 (2003) doi: 10.1038/nsb1203-980

Cite OneDep:

Structure **25**, 536-545 (2017) doi: 10.1016/j.str.2017.01.004

wwPDB Members




wwPDB.org





5 EASY STEPS TO PDB DEPOSITION





5 Easy Steps to PDB Deposition with OneDep

1 Check Sample Sequence

Sequence should

- Contain all residues
- Include expression tags and disordered residues

Check against reference database with BLAST

- Proteins ⇒
 UniProt BLAST
- Nucleic acids ⇒
 NCBI BLAST
- Any mismatch should correspond to mutation, variant or expression tag in your sequence.

Example of a UniProt alignment with highlighted mismatch

Query

61 RWWANDGRTPGSR

RWW NDGRTPGSR

79 RWWCNDGRTPGSR

P00698 (UniProt Sequence)

2 Check Ligands

Are your ligands already in our Chemical Component Dictionary (CCD)?

Check at:

- Ligand Expo ligand-expo.rcsb.org
- PDBeChem pdbe.org/chem
- Chemie search pdbj.org/chemie-search

Is the ligand present in the CCD?

- The 3-character code for your ligand should match that in the dictionary
- No During deposition you can provide a SMILES or upload a 2D chemical diagram

Example of a chemical diagram that should be provided at deposition

3 Prepare Data

Generate coordinate file in PDBx/mmCIF format

- This can be output directly from some refinement programs (e.g., PHENIX, REFMAC)
- This ensures that the maximum metadata is captured from the file for deposition


pdb_extract can be used
to extract information.
Visit pdb-extract.wwpdb.org

- Creates single PDBx/ mmCIF file containing extra meta data (e.g., author info, data collection and refinement stats) using output log files
- Converts SF files to PDBx/mmCIF format for non-MTZ format files

4 Validate! Validate!

The wwPDB provides a standalone validation server for checking your structure Visit validate.wwpdb.org

- Both model and experimental data are validated
- Provides detailed reports (PDF and XML files)
- Attempt to fix major issues highlighted prior to deposition to minimize the requirement for replacement files. It will reduce the time needed to process your entry.



Deposit your Structure

Deposition to the OneDep system can be carried out at the wwPDB deposition pages. Visit deposit.wwpdb.org

Deposition is tailored to experimental type:

- X-ray, Neutron, NMR, EM, or any combination of these
- You must upload the relevant coordinate and experimental data files and include mandatory information

A session ID is provided to enable you to continue your deposition at a later date. A PDB ID will be provided upon completion of your deposition.

After Deposition

Entries are annotated, validated, and returned to authors for review

 Processed files and official validation report can be downloaded at communication page

Request entry release and provide citation publication information using OneDep's communication panel.

Coordinate release follows instructions given at deposition

- REL released immediately after processing
- **HPUB** hold until publication, up to 1 year
- HOLD hold up to 1 year

If the structure is not published within 1 year, you have the option to either release or withdraw the entry.

